



Portable Autoclave (P-Type)



An ISO 9001:2000 Company

It is a unique vessel that generates pressure of 15 lbs and temperature of 121°C. Sterilization of any media is guaranteed with rapid speed.

The vessel is made of high grade aluminum/stainless steel (grade SS 304) with deep drawing bottom & top lid. Unit is provided with very efficient one turn locking system with perfect sealing. A stop cock dead weight valve is provided as a safety device.

Size : 12" x 12" Aluminium/Stainless Steel

Size : 12" x 14" Stainless Steel